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

Diafatur® - the new prep-technique


Drehzahl / speed
30'-300000 Upm / rpm

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1. Removing contact points, mesial-distal (without danger of injury).
2. Determination of prep boundary.
3. Circular grinding without damaging the prep-stage.
4. Precise grinding of the prep-stage with the help of a precision guide-pin.

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Pin not diamond-coated, avoids injury of the gingival edge.

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
Diafatur® FGs are packed as 2 pieces.
F = ISO-Shape-No. · hl = head length in mm
d = diameter in 1/10 mm · P = price category

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In case of all-ceramic restorations the minimum requirements in terms of wall strength must be met for stability reasons. Wall strengths at the lateral surfaces should be a minimum of 0.8 - 1 mm, at the crown margin 0.5 - 0.8 mm. Moreover, restorations using zirconium oxide and aluminium oxide as framework require preparation

with an accentuated shoulder at an angle of 45 - 90°, preferably with round shoulders. Arbitrary procedures are never tolerated by modern all-ceramic restorations, be they press or layer ceramics or zirconia.

Diafutor® C for silicate ceramics (feldspar and glass ceramics, LS2) and Diafutor® Z for oxide ceramic systems (Al_2O_3 , ZrO_2) take account of the required minimum strengths including a safety margin by virtue of their design. This design facilitates simple contouring of a 90° round shoulder. The precision guide pin helps to achieve a controlled removal of tooth material at a defined penetration depth, at the same time the pin also protects against injuries to the gingival margin.

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For all-ceramic restorations with silicate ceramics, e.g. LS2.

F :	508	001	152	509	509
d :	014	018	020	020	020
KL/hl :	8	-	8	8	8
P :	3	1	9	9	9

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For all-ceramic zirconia restorations.

F :	508	001	152	509	509
d :	014	014	016	018	018
KL/hl :	8	-	8	8	8
P :	3	1	3	3	3

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Diafutor® FGs are packed as 2 pieces.
F = ISO-Shape-No. · hl = head length in mm
d = diameter in 1/10 mm · P = price category