

Accessories for VECTOR® screws

Learn more:



Serrated housing nut

For single tooth movement using VECTOR® piston screws during an orthodontic treatment with removable appliances.

// Serrated housing nut

2481.1
5 pieces



Mounting screw

For fixing the serrated housing nut during an orthodontic treatment using VECTOR® piston screws.

Note: Mounting screw is not classified as medical product and not suitable for use in the mouth, contains nickel.

// Mounting screw

2480.1
5 pieces



Reamer

For recutting the thread when fabricating a removable appliance incorporating VECTOR® piston screws.

// Reamer

2482.1
1 piece



Screwdriver F

For activation of VECTOR® piston screws.

The screwdriver is not classified as medical product and not suitable for use in the mouth.

// Screwdriver F

2483.1
1 piece



Key for expansion screws

For precise adjustment of VECTOR® expansion screws.

// Key for expansion screws

2101.2
100 pieces



Security key for expansion screws

For HF5, HF9 and HF15 telescopes (IST® devices) and all kinds of palatal expansion screws.

// Security key for expansion screws

2100.1
10 pieces

