

Accessories

Learn more:



CLEAR-BLOKKER®

Light-curing, high-transparent gel for building-up pressure moulded DURAN® splints, reinforced with DURASPLINT® LC. Gel is best suited for removal of air inclusions and insertion of canine guidance. Polymerisation in wavelength range 360–420 nm, for example with LC Mini Light or LC-6 Light Oven. Prior to processing, LC primer has to be applied to DURAN® splints.

// 4 x 1.5 g
incl. 2 cannula Ø 0.4 mm and 8 cannula Ø 0.9 mm

5261.1
1 piece

CE



SD Cyano Veneer Fast / SD Pipettes

Fast-curing bonding agent, fluid and crystal clear, suitable for fixation of components when fabricating TAP® or IST® splints. Can also be used for secure adhesion of screws in sleep medicine and orthodontics.

Separately available SD pipettes ensure precise dosing of SD Cyano Veneer Fast by micro tip.

// SD-Cyano Veneer Fast

3701.1
5 g

// SD-Pipettes

3703.1
25 pieces

CE



Activation key

Activation key for TAP®-T and TAP®-T Reverse protrusion adjustment device. Can also be used for IST® hexagon screws.

// Activation key

9035.1
1 piece

