

# Bite registrations and occlusion trainer

Learn more:



## Bite Fix®

Bite registration with reusable holder made of sterilizable plastic for the simple and precise registration of construction bites in anti-snoring therapy. The decentralized handle assures easy midline control and comes with fixing screw and mm-scale for control of protrusion adjustment.

// kit 5474.1  
incl. 1 holder and 1 bite fork 1 kit

// bite forks 5475.1  
10 pieces



## Suitable for model scanners

## Bite Fix® Scan

Bite registration with reusable holder made of sterilizable plastic for the simple and precise registration of construction bites in anti-snoring therapy. The decentralized handle assures easy midline control and comes with fixing screw and mm-scale for control of protrusion adjustment. Fits into every model scanner due to a predetermined breaking point on the bite fork.

// kit 5492.1  
incl. 1 holder and 1 bite fork 1 kit

// bite forks 5493.1  
10 pieces



## New

## Bite Fix® One

Bite Fix® bite registration as single-use product, for easy and precise registration of construction bites in the anti-snoring therapy. Non-reusable, not sterilisable.

// Bite Fix® One 5503.1  
1 kit



## GEORGE GAUGE™ bite registration set

For registration of occlusion in snoring therapy and orthodontics. The bite registration is pre-adjusted using 2 or 5 mm bite forks.



// GEORGE GAUGE™ bite registration set comes with 3 bite forks (1 for each size)	<b>5471.1</b> 1 set
// GEORGE GAUGE™ bite fork, 5 mm, white, large (L)	<b>5462.1</b> 10 pieces
// GEORGE GAUGE™ bite fork, 2 mm, grey, small (S)	<b>5472.1</b> 10 pieces
// GEORGE GAUGE™ bite fork, 5 mm, white, small (S)	<b>5473.1</b> 10 pieces

CE

## AM Aligner

Thermoplastic occlusion trainer for condylar reposition after wearing protrusion splints during the night. The AM Aligner is heated and adapted to the patient's upper jaw by the clinician.



// AM Aligner	<b>9072.1</b> 10 pieces
---------------	----------------------------

CE

## Occlusion trainer

Condylar reposition device for patients' training after wearing protrusion splints during the night.



// Occlusion trainer	<b>5470.1</b> 10 pieces
----------------------	----------------------------

CE