



COMPOSITE INSTRUMENTS TITANIUM COATED

STRUMENTI PER COMPOSITO RICOPERTI IN NITRURO DI TITANIO



499/1T

(1:1 scale / scala 1:1)




1,9 mm




499/1T

2,1 mm




499/2T

2,2 mm



499/3T

2,2 mm




499/4T

ø 1,9 mm




499/5T

1,7 mm



499/6T

1,2 mm



499/7T

2,0 mm




499/8T

2,0 mm




499/9T

1,2 mm




499/11T




499/12T




2,1 mm




2,4 mm




ø 1,3 mm



ø 1,3 mm



ø 1 mm



ø 1 mm



1,3 mm




2,5 mm



2,5 mm




1,1 mm



499/SET

Composed of: 499/1T, 499/2T, 499/4T, 499/5T, 499/7T (contained can be changed under request)



“ MEDESY IN SPOT

The titanium coating of instruments (in yellow) minimizes the adhesion of composite.

La copertura in titanio degli strumenti (parte in giallo) riduce l'attaccamento del materiale composito. ”

1675/3

Composite posterior kit
Kit per composito posteriore

Contains / Contiene:

479/29	Burnisher White
499/1T	Composite instr. Titanium
499/5T	Composite instr. Titanium
499/6T	Composite instr. Titanium
499/9T	Composite instr. Titanium
521/1	Filling instr. Romerowski
575/1	Ceramic instr.
575/3	Ceramic instr.
581/1	Modelling instr. P-Thomas
980/9	Gammatrix tray Doi



RESTORATIVE / CONSERVATIVA

1675/4

Composite anterior kit
Kit per composito anteriore

Contains / Contiene:

499/3T	Composite instr. Titanium
499/6T	Composite instr. Titanium
499/8T	Composite instr. Titanium
499/11T	Composite instr. Titanium
499/12T	Composite instr. Titanium
575/1	Ceramic instr.
575/3	Ceramic instr.
581/1	Modelling instr. P-Thomas
980/9	Gammatrix tray Doi

494/2

(1:1 scale / scala 1:1)



2 mm

2,5 mm

2,9 mm




1 mm

1,45 mm

2 mm

1 mm

1,45 mm



518/1
SS1

518/2
SS2

518/3
SS3


494/1
BS1

494/2
BS2


494/3
BS3

497/1
BS4


497/2
BS5




2 mm




2,5 mm




2,9 mm



1,7 mm



2 mm



2,5 mm



1,7 mm



2 mm

2 mm

1,2 mm


1,7 mm

1,2 mm

1,7 mm




497/3
BS6



484/1
BP1



484/2
BP2




485/1
BP3



485/2
BP4



473/1
Woodson N. 1




473/2
Woodson N. 2




473/3
Woodson N. 3




2,5 mm




1,8 mm




2,4 mm





1,6 mm



1,8 mm





HEIDEMANN SPATULAS

SPATOLE DI HEIDEMANN

