This product is a Wax Pot, a kind of the dental technician tools that melts the solid wax for dental usage into the state of liquid.

- · A certain amount of 4 different waxes can be put and melted without the change of the property of the substance.
- · It is possible to check and compare the current temperature and set temperature.
- · It is a product that user can make precise technical tools through the maximum reducing of the temperature deviation range.
- · Control of function is convenient and is designed ergonomically for the convenient working by the user.
- · Wax is melted fast and the attribute of wax is not changed while the state of liquid is maintained.
- · Required temperature can be set correctly due to the Digital Circuit equipped with CPU.
- · Celsius and Fahrenheit (°C/°F) can be conveniently selected.

본 제품은 고체상태의 치과기공용 왁스를 액체상태로 녹이는 치기공기기의 일종인 Wax Pot이다.

- · 온도 설정 후 물성의 변화 없이 일정량의 4가지 왁스를 녹일 수 있다.
- · 현재온도와 설정온도를 비교하면서 작업이 가능하다.
- · 온도편차의 폭을 최대한 줄여 정밀한 기공물 제작이 가능하다.
- · 기능조작이 편리하며 사용자가 작업이 편리하도록 인체공학적으로 설계되었다.
- · 왁스가 빠르게 녹고 액체상태를 유지하면서 왁스의 속성이 변하지 않는다.
- · CPU를 장착한 Digital회로로써 원하는 온도를 정확하게 설정할 수 있다.
- · 섭씨와 화씨($^{\circ}$ C/ $^{\circ}$ F)를 편리하게 선택할 수 있다.

Specifications		
Classification		Input Power Source
DS-400		AC 100V-240V 50/60Hz
*Power 120W	*FUSE 4A	*Weight(kg) 1.13
*Temperature : 50°C~120°C, 122°F~248°F		
*Dimension(mm): W192 / D114 / H70		

