Finishing and polishing



SD-Foil scissors A

SD foil scissors A are delicately made of a high-grade auto-clavable stainless steel alloy. For better handling, scissors are laterally curved. Micro serration prevents slipping.

// SD-Foil scissors A

3460.1 1 piece

C



SD-Foil scissors B

SD foil scissors B are robust and durable for cutting splints or inserting cuts into splints after pressurising.

// SD-Foil scissors B

3461.1



Finishing set

Set containing eight grinding instruments for cutting, trimming and finishing of soft-elastic and hard-elastic material. With color coding for easy identification.

// Finishing set

3378.1

Set comes with one each cutting bur, finishing bur, HM carbide bur, ^{1 piece} HM tricutter fine, BIOPERM trimmer, POLYFIX, OSAMU polisher, DIMO®PRO SLIM and dressing stone.



LC-Bur

Specifically designed for finishing cold cure resin or light cured adjustments on splints. Allows neat cutting without clogging and leaving any residue. LC bur A comes in conical form for specific finishing or occlusal grinding. LC bur B is pear-shaped, ensuring very smooth grinding surface patterns when it comes to larger surfaces and LC bur C comes in bud form for quick and uniform removal of roughness.



Finishing bur

Square cutting bur made of special soft alloy that won't damage clasp wires.

HM-Tricutter fine

Fine triangular bur for interdental finishing of splints.

HM-Carbide cutter

Triangular carbide cutting bur for quick, neat cutting of hard pressure moulding material.

Cutting bur

Spiral cutting bur for quick, neat cutting of hard pressure moulding material.





OSAMU Polisher

Finishing and polishing of soft-elastic and compound material.

POLYFIX

For trimming and final polishing of soft-elastic and compound material, brown. The discs should be shaped before use with a sharpening stone.



For grinding and finishing soft-elastic and sandwich material, blue.

OSAMU Polisher, grey maximum rotational speed: 15,000 r/min	3247.1 2 pieces
// POLYFIX, on mandrel maximum rotational speed: 25,000 r/min	3371.1 1 piece
// POLYFIX, discs maximum rotational speed: 25,000 r/min	3372.1 6 pieces
✓ BIOPERM® Trimmer, blue maximum rotational speed: 25,000 r/min	3226.1 2 pieces



$DIMO^{\circledR}$

Stable and flexible trimming wheels with prolonged processing time, three types of grains, incl. mandrel.

Application

Trimming and initial polishing of elastomer, sandwich material, mouth guard and orthodontic appliances.

// brown / coarse, incl. 1 mandrel maximum rotational speed: 10,000 r/min	3382.1 10 pieces	3382.2 50 pieces
// green / medium, incl. 1 mandrel maximum rotational speed: 10,000 r/min	3383.1 10 pieces	3383.2 50 pieces
maroon / fine, incl. 1 mandrel maximum rotational speed: 10,000 r/min	3384.1 10 pieces	3384.2 50 pieces
// mandrel	3380.1 10 pieces	



DIMO® PRO / DIMO® PRO SLIM

Trimming wheel made of special elastomer with fine fibrous structure, extremely long working time, incl. mandrel. The discs should be shaped before use with a sharpening stone.

Application DIMO®PRO

For finishing large-scale marginal areas and for final and high-gloss polishing of splints, positioners, mouth guard and orthodontic appliances. Provides highly-polished rims and uniform surfaces of hard/soft compound material.

Application DIMO®PRO SLIM

For finishing splints, specifically in the area of ligaments and interdentally.

DIMO®PRO, Ø 24 x 6 mm, incl. 1 mandrel maximum rotational speed: 10,000 r/min	3381.1 10 pieces
// DIMO®PRO SLIM, Ø 24 x 4 mm, incl. 1 mandrel maximum rotational speed: 10,000 r/min	3376.1 10 pieces



High-quality goat hair brush with micro fibre insert, for use in combination with SD polishing pastes A and B for polishing adjusted splints and orthodontic

SD-Polishing mop for high gloss

Extra soft mop structure for finishing and high-gloss polishing of built-up splints and orthodontic appliances.

SD-Felt wheel

For use in combination with SD polishing pastes when polishing built-up splints and orthodontic appliances.

SD-Polishing brush, with mandrel maximum rotational speed: 5,000 r/min	3712.1 5 pieces
SD-Polishing mop for high gloss, with mandrel maximum rotational speed: 2,500 r/min	3714.1 5 pieces
// SD-Felt wheel, with mandrel maximum rotational speed: 5,000 r/min	3713.1 5 pieces

SD-Polishing paste

Specifically designed for finishing (SD polishing paste A) and for high-gloss polishing (SD polishing paste B) of splints, adjusted with DURASPLINT® or DURASPLINT® LC. Can also be used for polishing orthodontic appliances.

// A (for finishing)	3710.1 60 g
B (for high-gloss polishing)	3711.1 100 g





