

IST® and OPM

Learn more:





IST® Kit HF 9

Intra-oral anti-snoring device, pressure-moulded. Connection between upper and lower jaw is achieved by continuously adjustable HF telescope. Lower jaw and tongue are advanced by protrusion.

// IST® Kit HF 9 5461.1
1 kit
incl. 2 telescopes, 4 DURAN® (2.0 mm), 4 ISOFOLAN®, DURASPLINT® resin, 4 fastening elements with O-ring, 4 SCHEU-Anchor, Alastik chain, silicone tube, security key, SD box for patient, occlusion trainer, 2 CETRON® cleaning powder bags of 15 g, instructions and patient information

€ 0044



Telescopes

Telescopes with limit stop, available in 3 sizes (HF 5, HF 9 and HF 15), continuously adjustable.

// HF 5	5449.1 1 pair
// HF 9	5454.1 1 pair
// HF 15	5444.1 1 pair

€ 0044

	Minimum length	Maximum length	Stroke	Length to be activated	$\left(4 \times \frac{1}{4}\right)$
HF 5	20 mm	28 mm	5 mm	3 mm	0.25 mm
HF 9	25 mm	39 mm	9 mm	5 mm	0.25 mm
HF 15	31 mm	51 mm	15 mm	5 mm	0.25 mm

Positioning aid for IST® (12° / 0°)

For simple and correct positioning of the IST® fastening elements.

// Positioning aid for IST® (12° / 0°) 5469.1
1 piece



Fastening elements for IST® appliance

With stable retention, come with O-rings.

// Fastening elements for IST® appliance
incl. silicone tube

5468.1
4 pairs

CE 0044



Fastening elements for IST® and Herbst appliances

With stable retention, come with hexagon screws.

// Fastening elements for IST® and Herbst appliances
incl. silicone tube

5457.1
4 pairs

CE 0044



Hexagon socket screws

Spare screws for fastening elements.

// Hexagon socket screws
for IST® appliance

5458.1
4 pieces

CE 0044



O-Rings

Replacement O-rings for IST® appliance, transparent.

// O-Rings
for IST® appliance

5466.1
10 pieces

CE 0044



OPM Kit

For simple fabrication of a pressure moulded OPM protrusion splint in snoring therapy. Activation is done in 1 mm steps by replacing the plastic bars on both sides.

// OPM Kit
incl. 6 pairs of plastic bars (25 - 30 mm long), 4 fastening
elements with O-rings, 1 bag CETRON® cleaning powder of 15 g,
occlusion trainer, silicone protection tube

5495.1
1 kit

CE

