

# SCALPEL HANDLES, BLADES REMOVER

MANICI PER BISTURI, PINZE TOGLI LAME



**3637/11**

mm 110

**Short handle**  
Manico corto




**3637/16**

mm 160

**Long handle**  
Manico lungo



**New scalpel handle adjustable in 7 positions and 180°**  
Nuovo manico per bisturi orientabile su 7 posizioni e 180°



**3639**

mm 125


**Razor blade breaker and holder**  
Pinza rompi e porta lame



**3636**

mm 160

**Scalpel blades remover**  
Pinza rimuovi lame



## EXAMPLE OF APPLICATION


- 1 Blade to be seized here;
- 2 Then lift it here;
- 3 Finally pull it off in forward direction.

## MODO DI IMPIEGO

- 1 La lama presa in questa posizione;
- 2 Inserire la punta fra il manico e la lama;
- 3 Sollevare e quindi spingere verso l'alto.

# SCALPEL HANDLES – SCALPEL BLADES

MANICI PER BISTURI – LAME STERILI



**3630**  
mm 165  
N. 7

**3631**  
mm 145

**3632**  
mm 145

**3629**  
mm 125  
Bad Parker  
Millimetered  
Millimetrata

**3633**  
mm 130  
N. 3

## STERIL SCALPEL BLADES N.3 LAME PER BISTURI STERILI N.3




**3635/10**




**3635/11**



**3635/12**



**3635/12D**



**3635/13**



**3635/15**



**3635/15C**

Blades in package of 100 pcs  
Lame in confezioni da 100 pz

## SCALPEL HANDLE AND BLADES FOR MICROSURGERY MANICO E LAME PER MICROCHIRURGIA



**3634**  
mm 135  
N. 4



**3638**  
mm 135  
Scalpel handle for micro blades  
Manico lame per microchirurgia



**3638/64**



**3638/65**



**3638/67**



**3638/69**



**3638/Ti**  
mm 135  
Titanium scalpel handle for micro blades  
Manico lame per microchirurgia in titanio



SURGICAL / CHIRURGIA