Accessories



CA® Splint bags

Resealable minigrip splint bags for passing on the CA® splints to the patient, can be labelled individually.



5283.1 100 pieces

Delivery exclusively to certified CA® CLEAR ALIGNER partners.



CA® Splint case for patients

Handy zipped case for storing and carrying CA® splints. Resealable CA® splint bags for passing on the CA® splints to the patient, can be labelled individually. Splint remover for easy removal of pressure moulded splints.

CA® Splint case for patients

5282.1

comes with CA® Travel Case, 3 CETRON® cleansing powder bags of 15 g, 3 CA® splint bags, splint remover

Delivery exclusively to certified CA® CLEAR ALIGNER partners.



CA® Travel Case

Retainer box for up to two CA® CLEAR ALIGNER splints. Specific shape to fit easily into purse or luggage. Detachable lid with additional space for expansion key, two ventilation holes.



5303.1

Delivery exclusively to certified CA® CLEAR ALIGNER partners.



CA® FLASH

For insulation and easy removal when fabricating pressure moulded CA® splints. CA® FLASH is applied to ISOFLAN® foil prior to pressure moulding.



100 ml pump spray bottle, without propellant gas

5308.1 100 ml



CA® Grinding set

Grinding set for finishing CA® CLEAR ALIGNER splints, comes with six different instruments and dressing stone.





CA® Plastic buttons

Transparent plastic buttons allow placing elastics in the treatment of class II and III.

// CA® Plastic buttons	5301.1 4 pieces	5301.2 40 pieces	
		(6	



CA® Cyano Veneer Fast / CA® Pipettes

Crystal-clear, fast-curing bonding agent for adhesion of CA® plastic buttons on CA® splints, can also be used for secure adhesion of screws for Herbst and IST® appliances. CA® pipettes for precise dosing of adhesive.

✓ CA® Cyano Veneer Fast	5302.1 5 g
✓ CA® Pipettes	5300.1 25 pieces
	C€



CA® Coloured pencils

Coloured pencils for marking the ${\rm CA}^{\rm @}$ measuring points on the hard plaster models, with red and blue lead.



