

CA[®] TIP Pliers

Learn more:



CA[®] TIP 1

Large spherical tip allows insertion of retention points in CA[®] splints, for increase of friction. Stainless steel, sterilisable.

Warranty: 5 years.

CA[®] TIP 1
hemisphere: Ø 1.45 mm

5281.1
1 piece

CE



CA[®] TIP 2

Small spherical tip allows forming of activation points in CA[®] splints. Stainless steel, sterilisable.

Warranty: 5 years.

CA[®] TIP 2
hemisphere: Ø 1.0 mm

5291.1
1 piece

CE



CA[®] TIP 3

Pliers with oval plateau, enable creation of activation surfaces in CA[®] splints. Stainless steel, sterilisable.

Warranty: 5 years.

CA[®] TIP 3
oval plateau: 3.6 x 1.5 mm

5292.1
1 piece

CE




CA® TIP 4

Pliers for torque control with curved plateau, enable creation of labial /lingual activation surfaces in CA® splints. Stainless steel, sterilisable.

Warranty: 5 years.

CA® TIP 4
curved plateau: 4.0 x 2.5 mm

5295.1
1 piece



New

CA® TIP 5

Pliers for aligners with angled pliers head for correction of splint edges in case of gingival pressure marks. Stainless steel, sterilizable.

Warranty: 5 years.

CA® TIP 5

5299.1
1 piece



New

CA® TIP 6

Pliers for aligners for removal of activation points and surfaces inserted by means of CA® TIP 2, 3 and 4. Also suitable for finishing and smoothing of the splint surface. Stainless steel, sterilizable.

Warranty: 5 years.

CA® TIP 6

5323.1
1 piece

