


## Präparationskit Veneers Veneer prep tool kit




200125

Bestell-Nr. Order number	ISO Nummer ISO number	Fig.Nr shape no.	Körnung grit	ø 1/10 mm ø 1/10 mm	KL mm hl mm
531518ZF	806 314	198	524	018	9
531528ZF	806 314	198	524	020	9
531516ZF	806 314	198	514	018	9
531526ZF	806 314	198	514	020	9
531568ZF	806 314	508	524	022	8
532828	806 314	257	524	023	5
530505	806 314	552	534	018	8
531566ZF	806 314	508	514	022	9
532826	806 314	257	514	023	5
532578ZF	806 314	141	524	014	9
532588ZF	806 314	141	524	016	9
532576ZF	806 314	141	514	014	9
532586ZF	806 314	141	514	016	9

Präparations-Set zum sicheren Anlegen von Lagern für Veneers. Der Tiefenmarkierer ermöglicht mit den tiefenkongruenten Instrumenten mit Führungsstift eine Präparationstiefe von mindestens 0,4 mm. Die Stirnradien sind nicht sphärisch sondern bilden ein Kugelsegment und ermöglichen Schultern wenn nötig bis zu 90° als Grundlage von ästhetisch hochwertigen Versorgungen.

*Tool-kit for guided preparations of veneers. The depth-markers in combination with the guide-pin ensure a prep-depth of at least 0.4 mm. The tip surfaces of the cylindrical and tapered instruments with non-spherical but round tops enable chamfer-preparations up to 90°. With these tools esthetical, high-quality restorations become reality.*

## Handyfix



Best. Nr. 54000  
Order No. 54000

Best. Nr. 54002  
Order No. 54002

Best. Nr. 54001  
Order No. 54001

Best. Nr. 54003  
Order No. 54003

### Handyfix-Bohrerstände, einzeln

- 77 Steckplätze für 37 FGs + 40 WST s
- thermisch belastbar bis 134°C
- 57 x 55 x 10 mm

### Handyfix trays, single

- 77 slots for 37 FGs + 40 CA s
- thermal stability up to 134°C
- 57 x 55 x 10 mm