



SYRINGES FOR INTRALIGAMENTAL ANAESTHESIA (4962 – 4963 – 4964)

The intraligament anaesthesia technique is used for the treatment of individual teeth by means of an injection between the alveolar crest and the tooth root. With this technique the needle is inserted through the gingival sulcus into the periodontal space along the mesial or distal surface of the tooth: the anaesthetic solution must be injected slowly keeping a constant and moderate pressure to facilitate the distribution of the liquid. Contraindicated in presence of deep periodontal pockets and gingival acute infection.

SIRINGHE PER ANESTESIA INTRALIGAMENTARE (4962 – 4963 – 4964)

L'anestesia intraligamentare è una tecnica secondo la quale i denti vengono anestetizzati singolarmente, per mezzo di una iniezione tra l'alveolo osseo e la radice del dente. La modalità di applicazione è la seguente: s'introduce l'ago nel solco gengivale lungo la superficie mesiale o distale del dente e lentamente si inietta la soluzione anestetica con una pressione costante (e non eccessiva) dando il tempo al liquido di diffondersi. È controindicata in presenza di tasche parodontali ed infezioni acute della gengiva marginale.



4962/1

(ml 1,8 EU)
(0,2 ml/click)

4962/2

(ml 1,8 AM)
(0,2 ml/click)

4963/1

(ml 1,8 EU)
(0,07 ml/click)

4963/2

(ml 2,2 AM)
(0,07 ml/click)



INOX SYRINGES

SIRINGHE INOX



4960/INOX

(ml 1,8 EU/AM)

Stainless Steel
Acciaio Inox



4956/INOX

(ml 1,8 EU/AM)

Stainless Steel
Acciaio Inox